



MURRAY MEDICAL SOLUTIONS

# Broncoflex<sup>®</sup>

Single use sterile video bronchoscope

Flexible video-bronchoscopes

**Made in France**

Diameters 5.5mm and 3.9mm



Uses :

- bronchial exploration:
  - bronchoalveolar lavage and protected specimen brush sampling,
  - remove airway secretions,
  - make emergency diagnosis,
- percutaneous tracheostomy,
- difficult and selective intubation,
- placing and checking of double lumen tube (**Broncoflex<sup>®</sup> S**).

Advantages :

- do away with disinfection constraints,
- available at any time,
- training assistance,
- save and share information,
- keep control of costs.





# MURRAY MEDICAL SOLUTIONS

## Technical characteristics:

### Broncoflex® S    Broncoflex® M

Field of direction	0°	0°
Field angle	120°	120°
Field depth	3 - 80 mm	3 - 80 mm
External diameter of the distal tip	3.9 mm	5.5 mm
External diameter of the insertion tube	3.7 mm	5.5 mm
Functional length	605 mm	605 mm
Internal diameter of the working channel	1.4 mm	2.1 mm
Top/bottom tip deflection angle	160°/160°	160°/130°
Tip deflection radius	7.5 mm	9 mm
Sterilisation	Ethylene oxide	Ethylene oxide

**Ergonomically designed and Adaptable,** for practitioner's benefit

**Performance and Safety,** for your patient's benefit

**Sterile** and single use, it eliminates any risk of cross-contamination.

**Available at any time** of the day.



**Developed and designed by** Axess Vision Technology  
The Broncoflex® are **easy to use** and powerful for **unsurpassed image quality.**



## Contact information:



MURRAY MEDICAL SOLUTIONS

PO Box 337, Rose Bay NSW 2029

t 02 9398 9950

e [info@murraymedicalsolutions.com.au](mailto:info@murraymedicalsolutions.com.au)

w [murraymedicalsolutions.com.au](http://murraymedicalsolutions.com.au)



**axessvision**  
TECHNOLOGY