

Broncoflex™

Single-use Bronchoscope



Easy BAL Clip™

TSC The Surgical Company

— ENDOVISION —

Specialization through Innovation

Easy BAL Clip™

EASY SAMPLING SOLUTION

Easy BAL Clip™ is a single-use sampling solution for broncho-alveolar lavage procedure. Attached to the Broncoflex, it can be connected to multiple sizes of specimen containers.



COMPATIBILITY WITH SPECIMEN CONTAINER

Easy BAL Clip™ allows connection to multiple sizes of specimen containers. No limit in the sample size or number of sampling.



SINGLE OPERATOR SAMPLING SOLUTION

A physician can solo perform sampling, and switch to another container as needed.



IMPROVED SUCTION LINE MANAGEMENT

Top handle placement keeps specimen container upright to prevent sample loss.

TECHNICAL SPECIFICATIONS

MATERIAL	ABS-PC
LENGTH	90mm
WIDTH	75mm
HEIGHT	30mm
WEIGHT	15.5g
PART NUMBER	41010000

The Surgical Company Endovision
Axess Vision Technology
Zone de la Liodière - 6 rue de la Flotière
37300 Joué-lès-Tours (France)
Tel.: + 33 (0)2 47 34 32 90

<https://tsc-group.com/endovision>
www.broncoflex.com

